

Microscope slides and cover glasses

– good price and quality

New!



Microscope slides of soda-lime-glass, in accordance with DIN ISO 8037/1, with CE-mark, manufactured in compliance with the Directive 98/79 EG (IVD), dimensions appr. 76 x 26 mm, thickness appr. 1 mm, precleaned, ready-to-use, double wrapping with cellophane foil, packed in boxes of 50 slides.

Description	Art. No:	Quantity per unit	LC	Price EUR per unit	Block qty. units	Block Price EUR per unit
Cut edges, non-frosted	7695001	50	1S6	1,70	20	1,45
Cut edges, frosted one end	7695002	50	1S6	2,30	20	2,00
Ground edges, corners 90°, frosted one end	7695003	50	1S6	2,90	20	2,50
Ground edges, corners 45°	7695004	50	1S6	2,65	20	2,30
Cut edges, colour frosted one end white	7695005	50	1S6	2,95	20	2,55
Cut edges, colour frosted one end blue	7695006	50	1S6	2,95	20	2,55
Cut edges, colour frosted one end rose	7695007	50	1S6	2,95	20	2,55
Cut edges, colour frosted one end yellow	7695008	50	1S6	2,95	20	2,55
Cut edges, colour frosted one end green	7695009	50	1S6	2,95	20	2,55
Cut edges, colour frosted one end orange	7695010	50	1S6	2,95	20	2,55
Cut edges, colour frosted one end purple	7695011	50	1S6	2,95	20	2,55
Ground edges, corners 45°, colour frosted one end white	7695012	50	1S6	4,00	20	3,40
Ground edges, corners 45°, colour frosted one end blue	7695013	50	1S6	4,00	20	3,40
Ground edges, corners 45°, colour frosted one end rose	7695014	50	1S6	4,00	20	3,40
Ground edges, corners 45°, colour frosted one end yellow	7695015	50	1S6	4,00	20	3,40
Ground edges, corners 45°, colour frosted one end green	7695016	50	1S6	4,00	20	3,40
Ground edges, corners 45°, colour frosted one end orange	7695017	50	1S6	4,00	20	3,40
Ground edges, corners 45°, colour frosted one end purple	7695018	50	1S6	4,00	20	3,40

Subject to technical modification. All prices are exclusive of VAT. Our general terms and conditions shall apply.

Th. Geyer Group: Renningen · Hamburg · Berlin · Lohmar · Roskilde (DK) · Höxter www.thgeyer.com

Adhesive microscope slides of soda-lime-glass, in accordance with DIN ISO 8037/1, with CE-mark, manufactured in compliance with the Directive 98/79 EG (IVD), dimensions appr. 76 x 26 mm, thickness appr. 1 mm, precleaned, ready-to-use, in foil-packed plastic boxes of 72 slides.



Description	Art. No:	Quantity per unit	LC	Price EUR per unit	Block qty. units	Block Price EUR per unit
Ground edges, corners 90°, colour frosted one end white, (+)positive charged	7695019	72	1S6	17,30	20	14,75
Ground edges, corners 45°, colour frosted one end white, (+)positive charged	7695020	72	1S6	17,50	20	14,90
Ground edges, corners 90°, colour frosted one end white, silanised surface	7695021	72	1S6	18,10	20	15,40

Subject to technical modification. All prices are exclusive of VAT. Our general terms and conditions shall apply.

Cover glasses made of borosilicate glass 3.3, in accordance with DIN ISO 8255/1, manufactured in compliance with the Directive 98/79 EG (IVD), plastic hinged lid boxes, shrunk in aluminium protective foil.



Description	Art. No:	Quantity per unit	LC	Price EUR per unit	Block qty. units	Block Price EUR per unit
Thickness no. 1 (0,13 - 0,16 mm), size 15 x 15 mm	7695022	200	1S6	3,00	10	2,55
Thickness no. 1 (0,13 - 0,16 mm), size 18 x 18 mm	7695023	200	1S6	2,60	10	2,25
Thickness no. 1 (0,13 - 0,16 mm), size 20 x 20 mm	7695024	200	1S6	3,35	10	2,85
Thickness no. 1 (0,13 - 0,16 mm), size 22 x 22 mm	7695025	200	1S6	3,65	10	3,15
Thickness no. 1 (0,13 - 0,16 mm), size 24 x 24 mm	7695026	200	1S6	4,40	10	3,75
Thickness no. 1 (0,13 - 0,16 mm), size 21 x 26 mm	7695027	200	1S6	4,40	10	3,75
Thickness no. 1 (0,13 - 0,16 mm), size 24 x 32 mm	7695028	100	1S6	3,10	10	2,65
Thickness no. 1 (0,13 - 0,16 mm), size 24 x 40 mm	7695029	100	1S6	3,60	10	3,10
Thickness no. 1 (0,13 - 0,16 mm), size 24 x 50 mm	7695030	100	1S6	4,10	10	3,50
Thickness no. 1 (0,13 - 0,16 mm), size 24 x 60 mm	7695031	100	1S6	4,65	10	4,00
Thickness no. 1.5 (0,16 - 0,19 mm), for Shandon Automatic Coverslipper, size 24 x 50 mm	7695032	100	1S6	4,90	10	4,20

Subject to technical modification. All prices are exclusive of VAT. Our general terms and conditions shall apply.

Th. Geyer

Germany

Th. Geyer GmbH & Co. KG

Dornierstraße 4
D-71272 Renningen

Phone +49 7159 1637-0
Fax +49 7159 18417
E-Mail renningen@thgeyer.de

Th. Geyer Hamburg GmbH & Co. KG

Oehleckerring 14
D-22419 Hamburg

Phone +49 40 533214-0
Fax +49 40 5321004
E-Mail hamburg@thgeyer.de

Th. Geyer - Vertriebsbüro West -

Kirchstraße 1
D-53797 Lohmar

Phone +49 2246 10098-0
Fax +49 2246 10098-29
E-Mail lohmar@thgeyer.de

Th. Geyer Berlin GmbH

Dillenburg Straße 53
D-14199 Berlin

Phone +49 30 391100-1
Fax +49 30 3932069
E-Mail berlin@thgeyer.de

Denmark

Th. Geyer Danmark ApS

Industrivej 6A
DK-4000 Roskilde

Phone +45 4630 00 30
Fax +45 4630 00 69
E-Mail sales@thgeyer.dk